


TOPI-CLICK® MICRO™

Micro Quick-Fill Technology™
Makes Filling Fast & Simple



Note: Graduated Lines are for filling only and not for patients. **Top** black line of the plunger indicates **Fill Volume**. The embossed **Refill Badge** helps patients reorder medication in a timely manner.



Step 1: Filling
Syringe Tip-to-Tip Filling-
Screw Luer syringe to the Luer lock tip of the Micro™. **Requires no extra Luer Lock adaptor.**



EMP Jar Quick-fill Adaptor™ Filling- Screw on Topi-CLICK® Micro™ EMP Jar Quick-fill Adaptor™. Fill topical medications directly from mixing jar with this added Micro™ accessory.

Select one of 3 applicators



Note: Soft Angled Applicator™ placement requires tip to be aligned where lower slant of tip is on the same side as **Refill Badge of the Micro™ tube.**

Step 2: Applicator Tip Attachment -
Select Applicator Tip and press on by pushing straight down (important) against a hard surface.



Step 3: Priming the Micro Advance the actuator to engage plunger by manually turning/CLICKing base. The Actuator Bar should seat completely to the plunger before priming. **Actuator can also be advanced by using a mini cordless screwdriver to speed up the process (200 RPM Max). A Micro Driver Socket Bit™ can be purchased for faster actuator advancement with a cordless mini screwdriver.** (See Advancing Actuator Instructions.)




Step 4: Cap Placement
Screw on cap securely to protect tip and prevent evaporation.



Step 5: Label Application
Apply above the base opposite the "REFILL Badge" indicator.

Step 3a - Priming



Base

Step 3a: Priming
Prime the Micro™ by turning/CLICKing the base until contents are dispensed, then wipe off. Be sure to align Micro™ base with the triad tube.
Approximate Priming: 0.15 mL or 3 CLICKs.

Dispense/CLICK: 0.05 mL **Approx. Residual: 0.15 mL**
Dispenser Volume: Up to 9.75 mL **Approx. Priming: 0.15 mL or 3 CLICKs**
Total Delivery Volume: Up to 9.45 mL



TOPI-CLICK[®] MICRO[™] Micro Quick-Fill Technology[™] Makes Filling Fast & Simple

Instructions to Advance the Micro[™] Actuator with a Mini Cordless Screwdriver



Advancing the Actuator with a Mini Cordless Screwdriver

Actuator can be advanced by using a mini cordless screwdriver without a reversible trigger (200 RPM Max). Accidentally turning the base backwards could cause the Micro to slip in your hand.

The BLACK & DECKER[™] 3.6 Volt Pink Rechargeable Electric Screwdriver is recommended.

A **Micro Driver Socket Bit[™]** can also be purchased for faster actuator advancement with a cordless screwdriver. (See info below.)

Topi-CLICK[®] Micro Driver Socket Bit[™] - SKU#: TC0910

The triad-shaped Micro driver is designed to perfectly conform to the Micro Base. It should be used with a 1/4 inch hex bit power screw driver. We recommend power drivers with a speed of 200 RPM or less to reduce the likelihood of damage to the Micro.

Step 1: Insert the **Micro Driver** into the screwdriver socket until it is properly seated. Depending on the manufacturer and model, there is usually a light click that lets you know when it is properly seated in the tool.

Step 2: Rotate the **Driver** in a clockwise direction to advance the **Micro Actuator** up the tube.



Notes:

We recommend stopping the Actuator Bar 1 mL short of the target graduation (currently only for priming 6 mL or smaller prescriptions.)

Step 3: Run up the **Actuator** using the **Actuator Bar** as a visual indicator. The **Actuator Bar** should seat completely to the plunger.

