

Dual-Purpose Silicone Pad

Convenient Solutions

Tap-n-CLICK™ silicone pads were designed with the patient and the pharmacy in mind.

For the patient, these innovative pads allow for “one-handed” CLICK-Dosing™ of medication.

For the pharmacy, these pads provide a sound dampening surface that assists with the tapping out of air bubbles during filling of medication.



Tap-n-CLICK™ Pads

Simple, Smart and Convenient Topical Dosing Accessories

Tap-It



CLICK-It



For Pharmacies – When Filling

1. Fill the Topi-CLICK with medication to a selected volume.
2. Tap the Topi-CLICK on the pad to remove air bubbles.
3. Repeat steps one and two as desired.

For Patients – When Medicating

1. Remove cap from the Topi-CLICK.
2. Place Topi-CLICK on the Pad.
3. Turn the Super-CLICKer™ to left to CLICK.
4. Apply Dose once desired medication is dispensed.

Ordering Information

SKU	Description	MSRP
Tap-n-CLICK™		

HCTNT

Tap-CLICK Pad™ Multi-Use Silicone Pad, 5/Pack

\$25.00