### Topi-CLICK® 140

Pain Management / HRT / Scarring / Wound Care



#### The Next Generation in Topical Dosing Technology

# CLICK-Assist™ Feedback

The CLICK-Assist - Patient Feedback System is the next step in automated CLICK-dosing technology aimed to improve patient compliance and therapy outcomes.

A proprietary Super-CLICKer $^{\text{\tiny M}}$  mechanism is at the heart of this  $2^{\text{nd}}$  generation dosing system. This new technology provides stronger audible CLICKs, heavier tactile dosing sensations and improved visual indicators. As a result, the patient will enjoy a more responsive dosing experience.

# Ergonomic Design



A friction-free domed applicator pad with a 4Port™ dispensing design provides a comfortable experience when applying larger doses of medication (0.50mL/CLICK) directly to specific therapy zones. All Topi-CLICK® patient friendly devices reduce the risk the cross-contamination of medication to loved ones and pets since no medication comes in contact with the hands during application. Elegantly designed and easy to operate, Topi-CLICK® applicators rest comfortably in the hand and provide an repeatable and accurate dose each time.

Topi-CLICK 140 dispensers continue to raise the bar in topical dosing technology.







The Topi-CLICK® 140 comes wrapped in a textured polypropylene skin with polished visual markers. The [REFILL-Badge™] is now twice the size for even better refill awareness and guidance. A polished indicator window provides a full length view of the medication level for easy volume identification.

All new next generation Topi-CLICK® dispensers come with a recycle code mold on the base for easy recycling information.



### Proven dosing technology – millions and millions used yearly and growing.

TC35 vs. TC140

The all new Topi-CLICK® 140 dispensers are 4X larger than the original Topi-CLICK® 35. With a larger volume dose of 0.50mL/CLICK average, they are perfect for Pain Management, HRT, Scar Care and Wound Care.

The Difference	Topi-CLICK® 35	Topi-CLICK® 140
Volume	35 mL	140 mL
Dosage/CLICK	(0.25mL/CLICK ave.)	(0.50mL/CLICK ave.)
Application Ports	1, 3	4



Topi-CLICK® 35



Topi-CLICK® 140



## 5 Ordering Information



Order Online at www.DoseLogix.com

SKU	Description	mL / Click	Qty.	MSRP
TC140BLUE	Topi-CLICK® 140, Blue, 140 mL, (0.5 0mL/CLICK)	0.50 mL/CLICK	10/Pack	\$77.70
TC140PINK	Topi-CLICK® 140, Pink, 140 mL, (0.50 mL/CLICK	0.50 mL/CLICK	10/Pack	\$77.70
TC140BLACK	Topi-CLICK® 140, Black, 140 mL, (0.50 mL/CLICK	0.50 mL/CLICK	10/Pack	\$77.70
TC140WHITE	Topi-CLICK® 140,White, 14 0mL, (0.50 mL/CLICK)	0.50 mL/CLICK	10/Pack	\$77.70
TC140GOLD	Topi-CLICK® 140, Gold, 14 0mL, (0.50 mL/CLICK)	0.50 mL/CLICK	10/Pack	\$77.70

Minimum Order Quantities Apply

### The Next Generation in Topical Dosing Technology

